

EXAM, STERILE, SURGICAL HAND PROTECTION

P R O D U C T C A T A L O G

Our Gloves are Disposable. The People they Protect are Not.®



Sempermed is a globally integrated manufacturer that provides focused sales, marketing, and logistical support to those in need of hand protection.

www.SempermedUSA.com • Phone 1.800.366.9545 • Outside the U.S. 727.787.7250



GLOVE the PLANET
www.GlovethePlanet.com

sempermed®

YOUR GLOBAL GLOVE EXPERT™

With nearly 100 years of manufacturing experience, Sempermed has long been a pioneer within the disposable glove industry. From our humble beginnings to ranking among the world's largest medical and industrial glove manufacturers, we have always remained committed to supplying quality hand protection products.

As a leading disposable glove supplier, we are able to leverage the latest manufacturing technology in quality and innovation. We continuously develop and test products designed to better serve the needs of our customers to ensure Sempermed's hand protection products achieve optimal safety, reliability and quality.

Our Gloves are Disposable The People they Protect are Not.®



Best Touch®

LATEX with Aloe & Vitamin E

LATEX EXAMINATION GLOVES



Order #	Size	Length	Palm Width
BTLA101	XS	220mm	70mm +/- 10mm
BTLA102	S	220mm	80mm +/- 10mm
BTLA103	M	230mm	95mm +/- 10mm
BTLA104	L	230mm	111mm +/- 10mm
BTLA105	XL	230mm	115mm +/- 10mm

*Measurements represent most recent available information.



HOW WE GLOVE THE PLANET WITH BEST TOUCH® LATEX

- ✓ Inner boxes made with recycled material
- ✓ Outer cartons made with recycled material



Not applicable
where prohibited.

Best Touch® Latex with Aloe and Vitamin E provides superior protection and comfort with the added benefit of aloe and vitamin E. This exclusive formula allows for an enhanced inner coating which helps prevent your skin from becoming damaged or irritated. Additionally, the vitamin E acts as an antioxidant providing extra protection for your hands.

GLOVE COUNT

100 Gloves/Box (90 Gloves/Box XL)
10 Boxes/Case

All Sempermed gloves are single use only

Caution: This product contains natural rubber latex which may cause allergic reactions. Safe use of this product by or on latex sensitized individuals has not been established.

Our Gloves are Disposable. The People they Protect are Not.®



StarMed®

LATEX EXAMINATION GLOVES



Order #	Size	Length	Palm Width
SM101	XS	220mm	70mm +/- 10mm
SM102	S	220mm	80mm +/- 10mm
SM103	M	230mm	95mm +/- 10mm
SM104	L	230mm	111mm +/- 10mm
SM105	XL	230mm	115mm +/- 10mm

*Measurements represent most recent available information.



HOW WE GLOVE THE PLANET WITH STARMED® LATEX

- ✓ Inner boxes made with recycled material
- ✓ Outer cartons made with recycled material



Not applicable where prohibited.

StarMed® Latex powder-free exam gloves are designed to provide the perfect combination of protection, tactile sensitivity and value. Our advanced manufacturing technology and online chlorination process creates a glove that is strong, flexible and goes on with ease. StarMed® Latex provides maximum protection and comfort at a great value.

GLOVE COUNT

100 Gloves/Box

10 Boxes/Case

All Sempermed gloves are single use only

Caution: This product contains natural rubber latex which may cause allergic reactions. Safe use of this product by or on latex sensitized individuals has not been established.

SemperCARE®

LATEX EXAMINATION GLOVES



Order #	Size	Length	Palm Width
SCLT101	XS	255mm	76mm +/- 10mm
SCLT102	S	255mm	84mm +/- 10mm
SCLT103	M	255mm	94mm +/- 10mm
SCLT104	L	255mm	105mm +/- 10mm
SCLT105	XL	255mm	113mm +/- 10mm

*Measurements represent most recent available information.



HOW WE GLOVE THE PLANET WITH SEMPERCARE® LATEX

- ✓ Inner boxes made with recycled material
- ✓ Outer cartons made with recycled material



Not applicable where prohibited.

SemperCARE® Latex is the second generation of our popular co-polymer latex glove, boasting enhanced elasticity, tactile sensitivity and sensation. This glove has a slick inner coating for ease of donning and added comfort. At Sempermed®, we care about your comfort and protection.

GLOVE COUNT

100 Gloves/Box

10 Boxes/Case

All Sempermed gloves are single use only

Caution: This product contains natural rubber latex which may cause allergic reactions. Safe use of this product by or on latex sensitized individuals has not been established.

Our Gloves are Disposable. The People they Protect are Not.®



SEMPERFORCE®

BLACK NITRILE EXAMINATION GLOVES
THAT BLACK GLOVE® | QUALITY AND CONSISTENCY FIRST.®



Order #	Size	Length	Palm Width
BKNF102	S	230mm	85mm +/- 5mm
BKNF103	M	230mm	95mm +/- 5mm
BKNF104	L	230mm	108mm +/- 5mm
BKNF105	XL	230mm	113mm +/- 5mm
BKNF106	XXL	230mm	118mm +/- 5mm

*Measurements represent most recent available information.



HOW WE GLOVE THE PLANET WITH SEMPERFORCE®

- Inner boxes made with recycled material
- Outer cartons made with recycled material

SemperForce® Black Nitrile examination gloves by Sempermed are designed for demanding work environments where your safety and comfort are a must. A textured surface promises secure gripping capabilities. SemperForce® Black Nitrile gloves are ideal for automotive, industrial, tattoo, and cosmetology applications.

GLOVE COUNT

100 Gloves/Box (90 Gloves/Box XXL)
10 Boxes/Case

All Sempermed gloves are single use only

Your purchase of this product helps us proudly support:



SEMPERSHIELD®

EXTENDED CUFF

NITRILE EXAMINATION GLOVES



Order #	Size	Length	Palm Width
SSNF102	S	280mm	89mm +/- 10mm
SSNF103	M	280mm	99mm +/- 10mm
SSNF104	L	280mm	119mm +/- 10mm
SSNF105	XL	280mm	129mm +/- 10mm
SSNF106	XXL	280mm	135mm +/- 5mm



HOW WE GLOVE THE PLANET WITH SEMPERSHIELD® NITRILE

- Inner boxes made with recycled material
- Outer cartons made with recycled material

*Measurements represent most recent available information.

With extra thickness and an extended cuff, SemperShield® Nitrile powder-free exam gloves are designed for high-risk applications where your safety and protection are a must. A textured surface promises a secure grip without compromising your sensitivity. SemperShield® Nitrile offers all the benefits of a nitrile glove plus it's safe for use with chemotherapy drugs.

GLOVE COUNT

50 Gloves/Box
10 Boxes/Case

Tested for use with chemotherapy drugs.

Meets USP 800 requirements for PPE compliance for chemotherapy.

All Sempermed gloves are single use only



Our Gloves are Disposable. The People they Protect are Not.®



SEMPERSHIELD®

NITRILE EXAMINATION GLOVES



Order #	Size	Length	Palm Width
SSSC102	S	240mm	84mm +/- 3mm
SSSC103	M	240mm	94mm +/- 3mm
SSSC104	L	240mm	105mm +/- 3mm
SSSC105	XL	240mm	113mm +/- 3mm
SSSC106	XXL	240mm	123mm +/- 3mm

*Measurements represent most recent available information.



HOW WE GLOVE THE PLANET WITH SEMPERSHIELD® NITRILE

- ✓ Inner boxes made with recycled material
- ✓ Outer cartons made with recycled material

SemperShield® Nitrile standard cuff exam gloves are designed for applications where your safety and protection are a must. Extra thickness and a fully textured surface promises a secure grip in dry and wet environments. SemperShield® Nitrile offers all the benefits of a nitrile glove plus it's ESD tested and approved, which is safe for use with sensitive electronics or flammable liquids.

GLOVE COUNT

100 Gloves/Box (90 Gloves/Box XXL)
10 Boxes/Case

All Sempermed gloves are single use only



SEMPERSHIELD® CR

CHLOROPRENE EXAMINATION GLOVES



Order #	Size	Length	Palm Width
SSCR101	XS	240mm	76mm +/- 4mm
SSCR102	S	240mm	86mm +/- 4mm
SSCR103	M	240mm	95mm +/- 4mm
SSCR104	L	240mm	106mm +/- 4mm
SSCR105	XL	240mm	115

*Measurements represent most recent available information.



HOW WE GLOVE THE PLANET WITH SEMPERSHIELD® CR CHLOROPRENE

- ✓ Inner boxes made with recycled material
- ✓ Outer cartons made with recycled material

SemperShield® CR Chloroprene exam gloves are designed for applications where your comfort and protection are critical. Chloroprene offers the fit and feel similar to latex, with the strength and protection similar to nitrile. Textured fingertips give enhanced gripping capabilities in dry and wet environments. SemperShield® CR Chloroprene also features Smart Dispense Technology, Single Glove Delivery™ packaging method, which reduces wrinkles and allows for smooth dispensing.

GLOVE COUNT

100 Gloves/Box
10 Boxes/Case

All Sempermed gloves are single use only



Our Gloves are Disposable. The People they Protect are Not.®

SemperCARE® NITRILE

by Sempermed®
EXAMINATION GLOVES



Order #	Size	Length	Palm Width
NIPFT101	XS	245mm	76mm +/- 10mm
NIPFT102	S	245mm	84mm +/- 10mm
NIPFT103	M	245mm	94mm +/- 10mm
NIPFT104	L	245mm	105mm +/- 10mm
NIPFT105	XL	245mm	113mm +/- 10mm

*Measurements represent most recent available information.



HOW WE GLOVE THE PLANET WITH SEMPERCARE® NITRILE

- ✓ Inner boxes made with recycled material
- ✓ Outer cartons made with recycled material

GLOVE COUNT

100 Gloves/Box
10 Boxes/Case

Tested for use with chemotherapy drugs.

Meets USP 800 requirements for PPE compliance for chemotherapy.

All Sempermed gloves are single use only



SemperSure®

NITRILE EXAMINATION GLOVES



Order #	Size	Length	Palm Width
SUNF201	XS	244mm	76mm +/- 5mm
SUNF202	S	244mm	86mm +/- 5mm
SUNF203	M	244mm	96mm +/- 5mm
SUNF204	L	244mm	106mm +/- 5mm
SUNF205	XL	244mm	116mm +/- 5mm

*Measurements represent most recent available information.



HOW WE GLOVE THE PLANET WITH SEMPERSURE®

- ✓ Inner boxes made with recycled material
- ✓ Outer cartons made with recycled material
- ✓ 200 gloves (180 XL) per box instead of 100 (90 XL)

SemperSure® by Sempermed nitrile examination gloves are designed for those with skin sensitivities to regular nitrile gloves. SemperSure® accelerator-free gloves offer all the benefits of a nitrile glove, minus the sulfur-based chemical accelerators that cause Type I (latex) and Type IV (chemical) allergies.

GLOVE COUNT

200 Gloves/Box (180 Gloves/Box XL)
10 Boxes/Case

Tested for use with chemotherapy drugs.

Meets USP 800 requirements for PPE compliance for chemotherapy.

All Sempermed gloves are single use only



Our Gloves are Disposable. The People they Protect are Not.®

Tender Touch®

NITRILE EXAMINATION GLOVES



Order #	Size	Length	Palm Width
TTNF201	XS	249mm	70mm +/- 10mm
TTNF202	S	249mm	80mm +/- 10mm
TTNF203	M	249mm	95mm +/- 10mm
TTNF204	L	249mm	110mm +/- 10mm
TTNF205	XL	249mm	120mm +/- 10mm

*Measurements represent most recent available information.



HOW WE GLOVE THE PLANET WITH TENDER TOUCH®

- Inner boxes made with recycled material
- Outer cartons made with recycled material
- 200 gloves (180 XL) per box instead of 100 (90 XL)

We've delivered the right combination of fit, feel, tactile sensitivity, strength and value in a latex-free exam glove. And better yet, we've delivered it in our Smart Dispense Technology, Single Glove Delivery® packaging. This innovative packaging method reduces wrinkles and allows for more gloves per box, making Tender Touch® a smart choice.

GLOVE COUNT

200 Gloves/Box (180 Gloves/Box XL)
10 Boxes/Case

Tested for use with chemotherapy drugs.

Meets USP 800 requirements for PPE compliance for chemotherapy.

All Sempermed gloves are single use only



StarMed® ROSE

with Aloe & Vitamin E
NITRILE EXAMINATION GLOVES



Order #	Size	Length	Palm Width
SMNR201	XS	220mm	70mm +/- 10mm
SMNR202	S	230mm	80mm +/- 10mm
SMNR203	M	230mm	95mm +/- 10mm
SMNR204	L	230mm	110mm +/- 10mm
SMNR205	XL	230mm	120mm +/- 10mm

*Measurements represent most recent available information.



HOW WE GLOVE THE PLANET WITH STARMED® ROSE

- Inner boxes made with recycled material
- Outer cartons made with recycled material
- 200 gloves (180 XL) per box instead of 100 (90 XL)

StarMed® ROSE Nitrile with Aloe & Vitamin E has what it takes to brighten any environment, from ER to exam room. An aloe and vitamin E coating soothes hands and helps prevent skin irritation, while StarMed® ROSE's formulation provides enhanced protection, comfort and tactile sensitivity. It also features Smart Dispense Technology, Single Glove Delivery® packaging method, which reduces wrinkles, waste and allows for more gloves per box.

GLOVE COUNT

200 Gloves/Box (180 Gloves/Box XL)
10 Boxes/Case

All Sempermed gloves are single use only



Your purchase of this product helps us proudly support:



Our Gloves are Disposable. The People they Protect are Not.®



StarMed® PLUS™

NITRILE EXAMINATION GLOVES



Order #	Size	Length	Palm Width
SMNP301	XS	240mm	75mm +/- 5mm
SMNP302	S	240mm	85mm +/- 5mm
SMNP303	M	240mm	95mm +/- 5mm
SMNP304	L	240mm	106mm +/- 5mm
SMNP305	XL	240mm	116mm +/- 5mm

*Measurements represent most recent available information.



HOW WE GLOVE THE PLANET WITH STARMED® PLUS™

- ✓ Inner boxes made with recycled material
- ✓ Outer cartons made with recycled material
- ✓ 300 gloves (250 XL) per box instead of 100 (90 XL)

StarMed® PLUS™ nitrile gloves are lighter than ever, strong as always. Their innovative design allows users to experience the benefits of unrestricted hand motion because they are so soft and flexible. The proprietary nitrile formula used during the manufacturing process conforms to your hand while improving upon the overall strength of the final product. Packed 300 gloves per box, (250 XL) these latex-free exam gloves are the best choice for operational efficiencies.

GLOVE COUNT

300 Gloves/Box (250 Gloves/Box XL)
10 Boxes/Case

Your purchase of this product helps us proudly support:



Tested for use with chemotherapy drugs.

Meets USP 800 requirements for PPE compliance for chemotherapy.

All Sempermed gloves are single use only



StarMed® ULTRA NITRILE

NITRILE EXAMINATION GLOVES



Order #	Size	Length	Palm Width
SMTN251	XS	220mm	70mm +/- 10mm
SMTN252	S	220mm	80mm +/- 10mm
SMTN253	M	230mm	95mm +/- 10mm
SMTN254	L	230mm	110mm +/- 10mm
SMTN255	XL	230mm	120mm +/- 10mm

*Measurements represent most recent available information.



HOW WE GLOVE THE PLANET WITH STARMED® ULTRA NITRILE

- ✓ Inner boxes made with recycled material
- ✓ Outer cartons made with recycled material
- ✓ 250 gloves (225 XL) per box instead of 100 (90 XL)

StarMed® Ultra Nitrile by Sempermed are thin by design and provide a comfortable fit and feel without sacrificing protection and peace of mind. It also features Smart Dispense Technology, Single Glove Delivery® packaging method, which reduces wrinkles and allows for more gloves per box. With all that it has to offer, StarMed® Ultra Nitrile is a popular choice for users in a variety of environments.

GLOVE COUNT

250 Gloves/Box (225 Gloves/Box XL)
10 Boxes/Case

Tested for use with chemotherapy drugs.

Meets USP 800 requirements for PPE compliance for chemotherapy.

All Sempermed gloves are single use only



Our Gloves are Disposable. The People they Protect are Not.®

sempermed®
YOUR GLOBAL GLOVE EXPERT™

Sempermed is a globally integrated manufacturer that provides focused sales, marketing, and logistical support to those in need of hand protection.

www.SempermedUSA.com • Phone 1.800.366.9545 • Outside the U.S. 727.787.7250



GLOVE the PLANET®
www.GlovethePlanet.com

StarMed® Select

NITRILE EXAMINATION GLOVES

Order #	Size	Length	Palm Width
SMNS101	XS	240mm	76mm +/- 4mm
SMNS102	S	240mm	86mm +/- 4mm
SMNS103	M	240mm	98mm +/- 4mm
SMNS104	L	240mm	107mm +/- 4mm
SMNS105	XL	240mm	115mm +/- 4mm



HOW WE GLOVE THE PLANET WITH STARMED® SELECT

- ✓ Inner boxes made with recycled material
- ✓ Outer cartons made with recycled material

*Measurements represent most recent available information.



StarMed® Select nitrile, is the optimal glove solution where comfort, strength, and tactile sensitivity go hand-in-hand. These latex-free exam gloves are packed 100 gloves per box, an economical choice for a variety of applications and environments where hand protection is critical. It also features Smart Dispense Technology, Single Glove Delivery™ packaging method, which reduces wrinkles and allows for smooth dispensing.

GLOVE COUNT

100 Gloves/Box

10 Boxes/Case

All Sempermed gloves are single use only



StarMed® BLACK

NITRILE EXAMINATION GLOVES

Order #	Size	Length	Palm Width
SMNB101	XS	240mm	76mm +/- 4mm
SMNB102	S	240mm	86mm +/- 4mm
SMNB103	M	240mm	98mm +/- 4mm
SMNB104	L	240mm	107mm +/- 4mm
SMNB105	XL	240mm	115mm +/- 4mm

*Measurements represent most recent available information.



HOW WE GLOVE THE PLANET WITH STARMED® BLACK

- ✓ Inner boxes made with recycled material
- ✓ Outer cartons made with recycled material



StarMed® Black nitrile offers a stellar combination of comfort, strength, and tactile sensitivity. The black color hides stains caused from blood, grease, dirt, dyes, and ink. These latex-free exam gloves are packed 100 gloves per box, an economical choice for a variety of applications and environments. StarMed® Black also features Smart Dispense Technology, Single Glove Delivery™ packaging method, which reduces wrinkles and waste. Experience smooth glove dispensing and easy donning that is out of this world.

GLOVE COUNT

100 Gloves/Box

10 Boxes/Case

All Sempermed gloves are single use only



Our Gloves are Disposable. The People they Protect are Not.®



SemperCARE®

STRETCH

VINYL EXAMINATION GLOVES



Strength. Durability. Sensitivity. Fit. The SemperCare® Stretch Vinyl glove is designed as the ultimate solution for those who are at risk of becoming sensitive to natural rubber latex. Its cream color and new resin are an upgrade from regular vinyl gloves, giving the SemperCare® Stretch Vinyl glove the look and feel of latex without the latex allergens. It offers exceptional strength and protection, yet maintains maximum tactile sensitivity.

GLOVE COUNT
100 Gloves/Box
10 Boxes/Case

All Sempermed gloves are single use only



*Measurements represent most recent available information.



HOW WE GLOVE THE PLANET WITH SEMPERCARE® STRETCH VINYL

- ✓ Inner boxes made with recycled material
- ✓ Outer cartons made with recycled material

SemperCARE®

Synthetic

VINYL EXAMINATION GLOVES



SemperCare® Synthetic brand examination gloves provide a solution for professionals and their patients who are at risk of becoming sensitive to natural rubber latex. They maximize tactile sensitivity and provide a contoured fit, reducing hand fatigue.

GLOVE COUNT
100 Gloves
10 Boxes/Case

All Sempermed gloves are single use only



*Measurements represent most recent available information.



HOW WE GLOVE THE PLANET WITH SEMPERCARE® SYNTHETIC

- ✓ Inner boxes made with recycled material
- ✓ Outer cartons made with recycled material

Our Gloves are Disposable. The People they Protect are Not.®

StarMed® STERILE Nitrile

STERILE POWDER-FREE NITRILE EXAMINATION GLOVES



EXAMINATION GLOVES

Order #	Size	Length	Palm Width
SMSN102	S (6-7)	240mm	86mm +/- 4mm
SMSN103	M (7-8)	240mm	98mm +/- 4mm
SMSN104	L (8-9)	240mm	107mm +/- 4mm

*Measurements represent most recent available information.

StarMed® Sterile Nitrile offers superb fit, feel, tactile sensitivity, and barrier protection. Designed for critical environments that require a sterile glove, but not necessarily a surgical glove, StarMed® Sterile Nitrile can be especially effective for wound care, trauma, and burn care. Sterile exam gloves are also appropriate for OB/GYN, CCU, and ER departments. All this – with the added benefit of gamma sterilization.

GLOVE COUNT

50 Pairs/Box
4 Boxes/Case

All Sempermed gloves are single use only

sempermed® syntegra cr

SURGICAL GLOVES



Order #	Size	Length	Palm Width
SCR600	6	265mm	76mm +/- 6mm
SCR650	6.5	265mm	83mm +/- 6mm
SCR700	7	265mm	89mm +/- 6mm
SCR750	7.5	265mm	95mm +/- 6mm
SCR800	8	265mm	102mm +/- 6mm
SCR850	8.5	265mm	108mm +/- 6mm

*Measurements represent most recent available information.

For the latex sensitive, Sempermed® Syntegra CR offers a similar fit and feel to latex. These quality surgical gloves are manufactured with polychloroprene, also known as neoprene, a trusted material with over 70 years of proven performance. Healthcare professionals have come to rely on the non-latex comfort and function of polychloroprene. Excellent tactile sensitivity, while the textured surface and beaded cuff enhance gripping and protection.

GLOVE COUNT

40 Pairs/Box
6 Boxes/Case

All Sempermed gloves are single use only

Our Gloves are Disposable. The People they Protect are Not.®



SURGICAL GLOVES

Order #	Size	Length	Palm Width
SPFP550	5.5	289mm	72mm +/- 6mm
SPFP600	6	289mm	78mm +/- 6mm
SPFP650	6.5	289mm	85mm +/- 6mm
SPFP700	7	289mm	91mm +/- 6mm
SPFP750	7.5	289mm	97mm +/- 6mm
SPFP800	8	289mm	104mm +/- 6mm
SPFP850	8.5	289mm	11mm +/- 6mm
SPFP900	9	289mm	116mm +/- 6mm

Caution: This product contains natural rubber latex which may cause allergic reactions. Safe use of this product by or on latex sensitized individuals has not been established.



Surgeons' hands are the most delicate instruments in the surgical suite. They need to be prepared. And protected. That's why we developed Sempermed® Supreme surgical gloves. These unique gloves are anatomically designed to reduce hand fatigue and increase comfort. They are powder-free to eliminate aerosolization and skin irritation. A special multi-layer copolymer composition provides optimal protection during surgery.

GLOVE COUNT

50 Pairs/Box
6 Boxes/Case

All Sempermed gloves are single use only

sempermed® syntegra IR

SURGICAL GLOVES

Order #	Size	Length	Palm Width
SIR550	5.5	270mm	73mm +/- 3mm
SIR600	6	270mm	79mm +/- 3mm
SIR650	6.5	270mm	85mm +/- 3mm
SIR700	7	280mm	91mm +/- 3mm
SIR750	7.5	280mm	97mm +/- 3mm
SIR800	8	280mm	105mm +/- 3mm
SIR850	8.5	285mm	111mm +/- 3mm
SIR850	9	285mm	112mm +/- 3mm

*Measurements represent most recent available information.



Sempermed® Syntegra IR is a latex-free surgical glove made of synthetic polyisoprene (PI). Free from latex proteins, this polyisoprene glove is manufactured with an innovative process that reduces the risk for sensitization. Sempermed® Syntegra IR can be used universally in all surgical applications, and is particularly recommended for emergency care, eye surgery, and pediatrics. The fully anatomical shape with curved fingers corresponds with the natural, relaxed hand position, preventing hand fatigue. The multi-layer coating means that this glove can be donned easily and quickly, even when double gloving.

GLOVE COUNT

40 Pairs/Box
6 Boxes/Case

All Sempermed gloves are single use only

Our Gloves are Disposable. The People they Protect are Not.®