



APPROVED

NIOSH N95
42CFR Part 84
TC-84A-5530

APPROVED

APPROVED

NIOSH N95
42CFR Part 84
TC-84A-7447

APPROVED

EFFECTIVE PROTECTION OF MEDICAL PERSONNEL

NIOSH N95

MS6115L / MS8225

Disposable Particulate Respirator

FACE MASKS

NIOSH N95

MS6115L Cup-type N95 masks

Disposable Particulate Respirator

PFE
≥95%
for
0.3 micron particles
(MPPS)



Contour Fit Design

Adjustable Nosepiece

Comfortable Nose Foam

Easy Breathing Speaking

NIOSH N95 APPROVED

APPROVED
NIOSH N95
42CFR Part 84
TC-84A-5530
APPROVED



Mask Type	Standards	Filtration Efficiency *
Respirator Mask (PFE Standard)		N95
	USA: NIOSH (42 CFR 84)	PFE ≥ 95% for 0.3-micron particles

* (PFE - Particle Filtration or BFE - Bacteria Filtration Efficiency)

Respiratory Mask of Class N95 is suitable for filtering of the most penetrating particles of 0.3-micron size that can be found in the air.

0.1 micron: Particle Filtration Efficiency standard (PFE).

0.3 micron: Used to represent the most-penetrating particle size (MPPS), which is the most difficult size particle to capture.

3.0 microns: Bacteria Filtration Efficiency standard (BFE) for masks per EN 14683.

FEATURES MS6115L Cup-type N95 masks

Name	MS6115L Cup-type N95
Colour	White
Item Number	6115L
Weight/unit (oz) (g)	±0.3oz 8.68g
Outer layer	PP Spunbond
Marking	As per NIOSH
Soft foam nosepiece	Black PU foam
Dimensions (in)	4.72*5.51*2.09 in (±0.20in)
Dimensions (mm)	120*140*53 mm (± 5mm)
Description	N95 Particulate Respirator
Material	3-plys, PP spunbond outer layer, PP meltbown high filtration layer, PP needle-punched nonwoven inner layer
Elastic head band	Latex-free elastic-width 7mm
Inner layer	PP needle-punched nonwoven
Filter media	PP meltbown high filtration layer
Adjustable nose clip	Aluminium
Characteristics	2 CFR standard. TC#84A-5530 95% filtration effective against particulate aerosols free of oils Aluminum nose-clip tex-Free elastic ultrasonic welded to cup

ORDERING INFORMATION

MS6115L Cup-type N95 masks

	Boxes	Cases
Materials	350 gsm white cardboard	Brown 5-ply single wall
Dimensions (in) (mm)	5.31*4.92*7.87 (in) 135*125*200 (mm)	25.59*11.42*16.53 (in) 650*290*420 (mm)
6115L	Gross Weight: (lb) 12.46lb/cs Net: 7.69lb/cs (kg) 5.65kg/cs Net: 3.49kg/cs	
Package 20 masks/polybag/box * 20 boxes = 400 masks / case		

NIOSH N95

MS8225 Folded N95 masks

Disposable Particulate Respirator



Soft smooth lining

Considerable fit design

Adjustable nosepiece

Comfortable nose foam

Easy breathing speaking

PFE
≥95%
for
0.3 micron particles
(MPPS)

APPROVED
NIOSH N95
42CFR Part 84
TC-84A-7447

APPROVED

Mask Type	Standards	Filtration Efficiency *
Respirator Mask (PFE Standard)		N95
	USA: NIOSH (42 CFR 84)	PFE ≥ 95% for 0.3-micron particles

* (PFE - Particle Filtration or BFE - Bacteria Filtration Efficiency)

Respiratory Mask of Class N95 is suitable for filtering of the most penetrating particles of 0.3-micron size that can be found in the air.

0.1 micron: Particle Filtration Efficiency standard (PFE).

0.3 micron: Used to represent the most-penetrating particle size (MPPS), which is the most difficult size particle to capture.

3.0 microns: Bacteria Filtration Efficiency standard (BFE) for masks per EN 14683.

FEATURES MS8225 Folded N95 masks

Name	MS8225 Folded N95
Colour	White
Item Number	8225
Weight/unit (oz) (g)	±0.34oz 9.6g
Outer layer	PP needle-punched nonwoven
Marking	As per NIOSH
Soft foam nosepiece	Black rectangular foam
Dimensions (in)	4.53*6.30 in (± 0.20in)
Dimensions (mm)	115*160 mm (± 5mm)
Description	N95 Particulate Respirator
Material	4-plys, PP spunbond outer layer, PP meltbown high filtration layer, PP needle-punched nonwoven inner layer
Elastic head band	Inner layer 19g, PP Spunbond layer
Inner layer	Spandex headstraps
Filter media	PP needle-punched nonwoven
Adjustable nose clip	PP meltbown high filtration layer
Characteristics	Iron wire + plastic nose clip
	Conforms to NIOSH 42 CFR standard.
	TC#84A-7447 95% filtration effective against particulate aerosols free of oils. Adjustable plastic nose clip, Latex-Free elastic ultrasonic welded to cup

ORDERING INFORMATION MS8225 Folded N95 masks

	Boxes	Cases
Materials	350 gsm white cardboard	Brown 5-ply single wall
Dimensions (in) (mm)	5.12*4.72*6.30 (in) 130*120*160 (mm)	24.41*11.02*13.38 (in) 620*280*340 (mm)
8225	Gross Weight: (lb) 11.90lb/cs Net: 8.15lb/cs (kg) 5.4kg/cs Net: 3.7kg/cs	
Package 20 masks/polybag/box * 20 boxes = 400 masks / case		